

Regucalcin

Catalog No: #AP87459



Package Size: #AP87459-1 50ug #AP87459-2 100ug #AP87459-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Regucalcin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RC,Gluconolactonase,GNL,Senescence marker protein 30,SMP-30,Smp30
Accession No.	Uniprot ID:Q64374
Uniprot	Q64374
GeneID	19733;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rgn

Note: This product is for in vitro research use only