

Serum amyloid A-4 protein

Catalog No: #AP87464



Package Size: #AP87464-1 50ug #AP87464-2 100ug #AP87464-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serum amyloid A-4 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Constitutively expressed serum amyloid A protein,C-SAA,CSAA
Accession No.	Uniprot ID:P35542
Uniprot	P35542
GeneID	6291;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SAA4

Note: This product is for in vitro research use only