

Solute carrier family 22 member 1

Catalog No: #AP87469



Package Size: #AP87469-1 50ug #AP87469-2 100ug #AP87469-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 22 member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Organic cation transporter 1,rOCT1,Oct1
Accession No.	Uniprot ID:Q63089
Uniprot	Q63089
GeneID	24904;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc22a1

Note: This product is for in vitro research use only