

Securin-2

Catalog No: #AP87470



Package Size: #AP87470-1 50ug #AP87470-2 100ug #AP87470-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Securin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pituitary tumor-transforming gene 2 protein
Accession No.	Uniprot ID:Q9NZH5
Uniprot	Q9NZH5
GeneID	10744;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PTTG2

Note: This product is for in vitro research use only