

Phosphoenolpyruvate carboxykinase, cytosolic [GTP]

Catalog No: #AP87472



Package Size: #AP87472-1 50ug #AP87472-2 100ug #AP87472-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phosphoenolpyruvate carboxykinase, cytosolic [GTP]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PEPCK-C,PEPCK1
Accession No.	Uniprot ID:P35558
Uniprot	P35558
GeneID	5105;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PCK1

Note: This product is for in vitro research use only