

Nuclear receptor ROR-gamma

Catalog No: #AP87475



Package Size: #AP87475-1 50ug #AP87475-2 100ug #AP87475-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor ROR-gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Nuclear receptor RZR-gamma,Nuclear receptor subfamily 1 group F member 3,RAR-related orphan receptor C,Retinoid-related orphan receptor-gamma,Thymus orphan receptor,TOR,Nr1f3,Rorg,Thor
Accession No.	Uniprot ID:P51450
Uniprot	P51450
GeneID	19885;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rorc

Note: This product is for in vitro research use only