

Small nuclear ribonucleoprotein E

Catalog No: #AP87476



Package Size: #AP87476-1 50ug #AP87476-2 100ug #AP87476-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Small nuclear ribonucleoprotein E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	snRNP-E, Sm protein E, Sm-E
Accession No.	Uniprot ID:A1XQR9
Uniprot	A1XQR9
GeneID	100038330;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SNRPE

Note: This product is for in vitro research use only