

Protein S100-A9

Catalog No: #AP87477



Package Size: #AP87477-1 50ug #AP87477-2 100ug #AP87477-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein S100-A9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Calgranulin-B,Leukocyte L1 complex heavy chain,Migration inhibitory factor-related protein 14,MRP-14,S100 calcium-binding protein A9,Cagb,Mrp14
Accession No.	Uniprot ID:P31725
Uniprot	P31725
GeneID	20202;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S100a9

Note: This product is for in vitro research use only