

Calcipressin-2

Catalog No: #AP87487



Package Size: #AP87487-1 50ug #AP87487-2 100ug #AP87487-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcipressin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Calcineurin inhibitory protein ZAKI-4,Down syndrome candidate region 1-like protein 1,Myocyte-enriched calcineurin-interacting protein 2,MCIP2,Regulator of calcineurin 2,Dscr111,Zaki4
Accession No.	Uniprot ID:Q9JHG2
Uniprot	Q9JHG2
GeneID	53901;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rcan2

Note: This product is for in vitro research use only