

Roundabout homolog 4

Catalog No: #AP87489



Package Size: #AP87489-1 50ug #AP87489-2 100ug #AP87489-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Roundabout homolog 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Magic roundabout, UNQ421/PRO3674
Accession No.	Uniprot ID:Q8WZ75
Uniprot	Q8WZ75
GeneID	54538;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ROBO4

Note: This product is for in vitro research use only