

Sigma non-opioid intracellular receptor 1

Catalog No: #AP87491



Package Size: #AP87491-1 50ug #AP87491-2 100ug #AP87491-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sigma non-opioid intracellular receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Aging-associated gene 8 protein,SR31747-binding protein,SR-BP,Sigma 1-type opioid receptor,SIG-1R,OPRS1,SRBP,AAG8
Accession No.	Uniprot ID:Q99720
Uniprot	Q99720
GeneID	10280;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SIGMAR1

Note: This product is for in vitro research use only