

## Putative adenosylhomocysteinase 2

Catalog No: #AP87496



Package Size: #AP87496-1 50ug #AP87496-2 100ug #AP87496-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Putative adenosylhomocysteinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AdoHcyase 2,IP3R-binding protein released with inositol 1,4,5-trisphosphate,S-adenosyl-L-homocysteine hydrolase 2,S-adenosylhomocysteine hydrolase-like protein 1,Irbit
Accession No.	Uniprot ID:Q80SW1
Uniprot	Q80SW1
GeneID	229709;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ahcyl1

Note: This product is for in vitro research use only