

Rho-associated protein kinase 1

Catalog No: #AP87497



Package Size: #AP87497-1 50ug #AP87497-2 100ug #AP87497-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rho-associated protein kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Liver regeneration-related protein LRRG199,Rho-associated, coiled-coil-containing protein kinase 1,Rho-associated, coiled-coil-containing protein kinase I,ROCK-I,p150 RhoA-binding kinase ROK beta,p160 ROCK-1,p160ROCK,Ac2-154
Accession No.	Uniprot ID:Q63644
Uniprot	Q63644
GeneID	81762;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rock1

Note: This product is for in vitro research use only