

Sphingosine 1-phosphate receptor 1

Catalog No: #AP87499



Package Size: #AP87499-1 50ug #AP87499-2 100ug #AP87499-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine 1-phosphate receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	S1P receptor 1,Endothelial differentiation G-protein coupled receptor 1,Lysophospholipid receptor B1,Sphingosine 1-phosphate receptor Edg-1,S1P receptor Edg-1,CD363,Edg1,Lpb1
Accession No.	Uniprot ID:O08530
Uniprot	O08530
GeneID	13609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S1pr1

Note: This product is for in vitro research use only