

Solute carrier family 22 member 1

Catalog No: #AP87501



Package Size: #AP87501-1 50ug #AP87501-2 100ug #AP87501-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 22 member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Organic cation transporter 1, OCT1
Accession No.	Uniprot ID:A7MBE0
Uniprot	A7MBE0
GeneID	519461;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC22A1

Note: This product is for in vitro research use only