

## Syndecan-2

Catalog No: #AP87502



Package Size: #AP87502-1 50ug #AP87502-2 100ug #AP87502-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Syndecan-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SYND2,Fibroglycan,Heparan sulfate proteoglycan core protein,HSPG,CD362,Hspg1,Synd2
Accession No.	Uniprot ID:P34900
Uniprot	P34900
GeneID	25615;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sdc2

Note: This product is for in vitro research use only