

## Phosphotriesterase-related protein

Catalog No: #AP87508



Package Size: #AP87508-1 50ug #AP87508-2 100ug #AP87508-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphotriesterase-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Parathion hydrolase-related protein,Resiniferotoxin-binding phosphotriesterase-related protein,Rpr-1
Accession No.	Uniprot ID:Q63530
Uniprot	Q63530
GeneID	63852;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pter

Note: This product is for in vitro research use only