

## Pleiotrophin

Catalog No: #AP87516



Package Size: #AP87516-1 50ug #AP87516-2 100ug #AP87516-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pleiotrophin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PTN,Heparin-binding brain mitogen,HBBM,Heparin-binding growth factor 8,HBGF-8,Heparin-binding growth-associated molecule,HB-GAM,Heparin-binding neutrophilic factor,HBNF,Osteoblast-specific factor 1,OSF-1
Accession No.	Uniprot ID:P63089
Uniprot	P63089
GeneID	19242;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ptn

Note: This product is for in vitro research use only