

Sphingosine 1-phosphate receptor 1

Catalog No: #AP87518



Package Size: #AP87518-1 50ug #AP87518-2 100ug #AP87518-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine 1-phosphate receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	S1P receptor 1, Endothelial differentiation G-protein coupled receptor 1, Sphingosine 1-phosphate receptor Edg-1, S1P receptor Edg-1, CD363, EDG1
Accession No.	Uniprot ID:Q5E9P3
Uniprot	Q5E9P3
GeneID	281135;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S1PR1

Note: This product is for in vitro research use only