

# RNA exonuclease 1 homolog

Catalog No: #AP87523



Package Size: #AP87523-1 50ug #AP87523-2 100ug #AP87523-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RNA exonuclease 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Transcription elongation factor B polypeptide 3-binding protein 1,Kiaa1138,Tceb3bp1
Accession No.	Uniprot ID:Q7TT28
Uniprot	Q7TT28
GeneID	66932;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Rexo1

Note: This product is for in vitro research use only