

Ribonuclease 7

Catalog No: #AP87524



Package Size: #AP87524-1 50ug #AP87524-2 100ug #AP87524-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribonuclease 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	RNase 7, Skin-derived antimicrobial protein 2, SAP-2, UNQ2516/PRO6006
Accession No.	Uniprot ID:Q9H1E1
Uniprot	Q9H1E1
GeneID	84659;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RNASE7

Note: This product is for in vitro research use only