

MAU2 chromatid cohesion factor homolog

Catalog No: #AP87526



Package Size: #AP87526-1 50ug #AP87526-2 100ug #AP87526-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAU2 chromatid cohesion factor homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAU-2,Cohesin loading complex subunit SCC4 homolog,KIAA0892,SCC4
Accession No.	Uniprot ID:Q9Y6X3
Uniprot	Q9Y6X3
GeneID	23383;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAU2

Note: This product is for in vitro research use only