

Kit ligand

Catalog No: #AP87527



Package Size: #AP87527-1 50ug #AP87527-2 100ug #AP87527-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kit ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Mast cell growth factor, MGF, Stem cell factor, SCF, c-Kit ligand, MGF
Accession No.	Uniprot ID:Q06220
Uniprot	Q06220
GeneID	403507;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KITLG

Note: This product is for in vitro research use only