

Sialic acid-binding Ig-like lectin 7

Catalog No: #AP87534



Package Size: #AP87534-1 50ug #AP87534-2 100ug #AP87534-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sialic acid-binding Ig-like lectin 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Siglec-7, Adhesion inhibitory receptor molecule 1, AIRM-1, CDw328, D-siglec, QA79 membrane protein, p75, CD328, AIRM1
Accession No.	Uniprot ID: Q9Y286
Uniprot	Q9Y286
GeneID	27036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SIGLEC7

Note: This product is for in vitro research use only