

DNA-directed RNA polymerase III subunit RPC9

Catalog No: #AP87535



Package Size: #AP87535-1 50ug #AP87535-2 100ug #AP87535-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA-directed RNA polymerase III subunit RPC9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	RNA polymerase III subunit C9, Calcitonin gene-related peptide-receptor component protein, CGRP-RCP
Accession No.	Uniprot ID:O75575
Uniprot	O75575
GeneID	27297;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CRCP

Note: This product is for in vitro research use only