

R-spondin-1

Catalog No: #AP87536



Package Size: #AP87536-1 50ug #AP87536-2 100ug #AP87536-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	R-spondin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cysteine-rich and single thrombospondin domain-containing protein 3,Cristin-3,Roof plate-specific spondin-1
Accession No.	Uniprot ID:Q9Z132
Uniprot	Q9Z132
GeneID	192199;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rspo1

Note: This product is for in vitro research use only