

Ras-related protein Rab-27B

Catalog No: #AP87540



Package Size: #AP87540-1 50ug #AP87540-2 100ug #AP87540-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ras-related protein Rab-27B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q99P74
Uniprot	Q99P74
GeneID	84590;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rab27b

Note: This product is for in vitro research use only