

cGMP-specific 3',5'-cyclic phosphodiesterase

Catalog No: #AP87542



Package Size: #AP87542-1 50ug #AP87542-2 100ug #AP87542-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	cGMP-specific 3',5'-cyclic phosphodiesterase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	cGMP-binding cGMP-specific phosphodiesterase,CGB-PDE,Pde5
Accession No.	Uniprot ID:Q8CG03
Uniprot	Q8CG03
GeneID	242202;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pde5a

Note: This product is for in vitro research use only