

Podoplanin

Catalog No: #AP87543



Package Size: #AP87543-1 50ug #AP87543-2 100ug #AP87543-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Podoplanin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Aggrus, Glycoprotein 38, Gp38, OTS-8, PA2.26 antigen, T1-alpha, T1A, Transmembrane glycoprotein E11, Gp38, Ots8 |
| Accession No. | Uniprot ID:Q62011 |
| Uniprot | Q62011 |
| GeneID | 14726; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: Pdpn

Note: This product is for in vitro research use only