

Protein SAAL1

Catalog No: #AP87544



Package Size: #AP87544-1 50ug #AP87544-2 100ug #AP87544-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein SAAL1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q96ER3
Uniprot	Q96ER3
GeneID	113174;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SAAL1

Note: This product is for in vitro research use only