

Syndecan-2

Catalog No: #AP87545



Package Size: #AP87545-1 50ug #AP87545-2 100ug #AP87545-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syndecan-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SYND2,Fibroglycan,Heparan sulfate proteoglycan core protein,HSPG,CD362,Hspg1,Synd2
Accession No.	Uniprot ID:P43407
Uniprot	P43407
GeneID	15529;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sdc2

Note: This product is for in vitro research use only