

# Scavenger receptor class B member 1

Catalog No: #AP87546

Package Size: #AP87546-1 50ug #AP87546-2 100ug #AP87546-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Scavenger receptor class B member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SRB1,CD36 and LIMPII analogous 1,CLA-1,CD36 antigen-like 1,Collagen type I receptor, thrombospondin receptor-like 1,SR-BI,CD36,CD36L1,CLA1
Accession No.	Uniprot ID:Q8WTV0
Uniprot	Q8WTV0
GeneID	949;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SCARB1

Note: This product is for in vitro research use only