

Scavenger receptor class B member 1

Catalog No: #AP87547



Package Size: #AP87547-1 50ug #AP87547-2 100ug #AP87547-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scavenger receptor class B member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	SRB1,High density lipoprotein receptor SR-BI,SR-BI
Accession No.	Uniprot ID:Q8SQC1
Uniprot	Q8SQC1
GeneID	397018;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SCARB1

Note: This product is for in vitro research use only