

Serpín H1

Catalog No: #AP87548



Package Size: #AP87548-1 50ug #AP87548-2 100ug #AP87548-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serpín H1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Collagen-binding protein, Colligin
Accession No.	Uniprot ID:Q2KJH6
Uniprot	Q2KJH6
GeneID	510850;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SERPINH1

Note: This product is for in vitro research use only