

Receptor-binding cancer antigen expressed on SiSo cells

Catalog No: #AP87549

Package Size: #AP87549-1 50ug #AP87549-2 100ug #AP87549-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-binding cancer antigen expressed on SiSo cells
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cancer-associated surface antigen RCAS1, Estrogen receptor-binding fragment-associated gene 9 protein, Rcas1
Accession No.	Uniprot ID:Q9D0V7
Uniprot	Q9D0V7
GeneID	55960;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ebag9

Note: This product is for in vitro research use only