

Endophilin-A2

Catalog No: #AP87551

Package Size: #AP87551-1 50ug #AP87551-2 100ug #AP87551-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endophilin-A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Endophilin-2, SH3 domain protein 2B, SH3 domain-containing GRB2-like protein 1
Accession No.	Uniprot ID:Q2KJA1
Uniprot	Q2KJA1
GeneID	512719;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SH3GL1

Note: This product is for in vitro research use only