

Semaphorin-3A

Catalog No: #AP87552



Package Size: #AP87552-1 50ug #AP87552-2 100ug #AP87552-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Semaphorin-3A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Semaphorin III,Sema III
Accession No.	Uniprot ID:Q63548
Uniprot	Q63548
GeneID	29751;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sema3a

Note: This product is for in vitro research use only