

Mothers against decapentaplegic homolog 1

Catalog No: #AP87553



Package Size: #AP87553-1 50ug #AP87553-2 100ug #AP87553-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAD homolog 1,Dwarfin-A,Dwf-A,Mothers-against-DPP-related 1,Mad-related protein 1,SMAD family member 1,SMAD 1,Madh1,Madr1
Accession No.	Uniprot ID:P70340
Uniprot	P70340
GeneID	17125;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Smad1

Note: This product is for in vitro research use only