

Mothers against decapentaplegic homolog 4

Catalog No: #AP87554



Package Size: #AP87554-1 50ug #AP87554-2 100ug #AP87554-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MAD homolog 4, Deletion target in pancreatic carcinoma 4 homolog, SMAD family member 4, SMAD 4, Dpc4, Madh4
Accession No.	Uniprot ID:P97471
Uniprot	P97471
GeneID	17128;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smad4

Note: This product is for in vitro research use only