

Serotonin N-acetyltransferase

Catalog No: #AP87555



Package Size: #AP87555-1 50ug #AP87555-2 100ug #AP87555-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serotonin N-acetyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Serotonin acetylase,Aralkylamine N-acetyltransferase,AA-NAT,Snat
Accession No.	Uniprot ID:O88816
Uniprot	O88816
GeneID	11298;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aanat

Note: This product is for in vitro research use only