

## Tyrosine-protein phosphatase non-receptor type 1

Catalog No: #AP87557



Package Size: #AP87557-1 50ug #AP87557-2 100ug #AP87557-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tyrosine-protein phosphatase non-receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein-tyrosine phosphatase 1B,PTP-1B,PTP1B
Accession No.	Uniprot ID:P18031
Uniprot	P18031
GeneID	5770;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PTPN1

Note: This product is for in vitro research use only