

Heterogeneous nuclear ribonucleoprotein A1

Catalog No: #AP87558



Package Size: #AP87558-1 50ug #AP87558-2 100ug #AP87558-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	hnRNP A1, Helix-destabilizing protein, HDP, Single-strand RNA-binding protein, hnRNP core protein A1, Hnrpa1
Accession No.	Uniprot ID:P04256
Uniprot	P04256
GeneID	29578;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hnrpa1

Note: This product is for in vitro research use only