

# Radical S-adenosyl methionine domain-containing protein 2

Catalog No: #AP87560

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP87560-1 50ug #AP87560-2 100ug #AP87560-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Radical S-adenosyl methionine domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Viperin, Virus inhibitory protein, endoplasmic reticulum-associated, interferon-inducible, Vig1
Accession No.	Uniprot ID:Q8CBB9
Uniprot	Q8CBB9
GeneID	58185;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Rsad2

Note: This product is for in vitro research use only