

Syntenin-1

Catalog No: #AP87563



Package Size: #AP87563-1 50ug #AP87563-2 100ug #AP87563-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syntenin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Syndecan-binding protein 1
Accession No.	Uniprot ID:Q9JI92
Uniprot	Q9JI92
GeneID	83841;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sdcbp

Note: This product is for in vitro research use only