

Semaphorin-6C

Catalog No: #AP87564



Package Size: #AP87564-1 50ug #AP87564-2 100ug #AP87564-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Semaphorin-6C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Semaphorin-Y, Sema Y, KIAA1869, SEMAY
Accession No.	Uniprot ID:Q9H3T2
Uniprot	Q9H3T2
GeneID	10500;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SEMA6C

Note: This product is for in vitro research use only