

Semaphorin-7A

Catalog No: #AP87565



Package Size: #AP87565-1 50ug #AP87565-2 100ug #AP87565-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Semaphorin-7A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Semaphorin-K1,Sema K1,Semaphorin-L,Sema L,CD108,Cd108,Semal,Semk1
Accession No.	Uniprot ID:Q9QUR8
Uniprot	Q9QUR8
GeneID	20361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sema7a

Note: This product is for in vitro research use only