

Sodium-dependent dopamine transporter

Catalog No: #AP87568



Package Size: #AP87568-1 50ug #AP87568-2 100ug #AP87568-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium-dependent dopamine transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DA transporter,Solute carrier family 6 member 3,Dat,Dat1
Accession No.	Uniprot ID:Q61327
Uniprot	Q61327
GeneID	13162;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc6a3

Note: This product is for in vitro research use only