

Rho-related GTP-binding protein RhoD

Catalog No: #AP87569



Package Size: #AP87569-1 50ug #AP87569-2 100ug #AP87569-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rho-related GTP-binding protein RhoD
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Rho-related protein HP1,RhoHP1,ARHD
Accession No.	Uniprot ID:O00212
Uniprot	O00212
GeneID	29984;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RHOD

Note: This product is for in vitro research use only