

# Septin-5

Catalog No: #AP87574



Package Size: #AP87574-1 50ug #AP87574-2 100ug #AP87574-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Septin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cell division control-related protein 1,CDCrel-1,Peanut-like protein 1,Pnutl1
Accession No.	Uniprot ID:Q9JJM9
Uniprot	Q9JJM9
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:43713

Note: This product is for in vitro research use only