

## Pulmonary surfactant-associated protein A1

Catalog No: #AP87589



Package Size: #AP87589-1 50ug #AP87589-2 100ug #AP87589-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pulmonary surfactant-associated protein A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PSP-A,35 kDa pulmonary surfactant-associated protein,Alveolar proteinosis protein,Collectin-4,COLEC4,PSAP,SFTP1,SFTPA,SFTPA1B
Accession No.	Uniprot ID:Q8IWL2
Uniprot	Q8IWL2
GeneID	653509;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SFTPA1

Note: This product is for in vitro research use only